



# README

**Objective:** Segment and extract arteries in the angiogram


**Processing pipeline:**




- Preprocessing:
  - Recreate the video for 21s. So, the number of frames will be 504
  - Crop each frame in order to obtain image with width=700 and height=600




- Image processing:
  - Apply *gaussian filter* to each frame (to smooth the image)




- Apply *frangi filter* and crop function (640,540) to each frame (to detect edges and remove noises)




- Apply *fillPoly* and *bitwise\_and* functions to remove noises from each frame



- Apply erosion with 5\*5 kernel to each frame (to further remove noises from each frame)




- Again, apply *gaussian* filter to each frame (to smooth the image)




- Video creation:
  - Integrate all the frames to create the video

#### Evaluation:

Original input



Final output



In final output, main components of the original image can be seen

Note: Example image is taken from,

[110.jpg](#)

**Run the code:**

- Import the shared `angiogram.mp4` to google colabs and run `angiogram_vedio_processing.ipynb` file

**References:**

- *Frangi filter*: [https://scikit-image.org/docs/0.14.x/auto\\_examples/filters/plot\\_frangi.html](https://scikit-image.org/docs/0.14.x/auto_examples/filters/plot_frangi.html)
- *Gaussian filter*: [https://opencv-python-tutroals.readthedocs.io/en/latest/py\\_tutorials/py\\_imgproc/py\\_filtering/py\\_filtering.html](https://opencv-python-tutroals.readthedocs.io/en/latest/py_tutorials/py_imgproc/py_filtering/py_filtering.html)